

GE 医疗中国

红色旋风
再次来袭

手动到自动…
蛋白纯化更简单



ÄKTA™ start

简单易学、简单易用的蛋白纯化系统，
把您从繁琐的手动纯化中解放出来

- 简单的纯化标签蛋白和抗体
- 获得实时监测
- 评估和分享你的结果





标签蛋白

通过触摸屏简单快速的直接控制纯化

缓冲液交换/脱盐

快速有效的组分分离且收率高——
使用UV280找到你的蛋白

抗体纯化

使用两个缓冲液入口方便进行结合和洗脱——使用样品入口阀简单的进行大体积样品的上样

非标签蛋白

离子交换和亲和层析（可以选择阶梯和线性洗脱）的方法模板——脱盐和分子排阻也是可以选的

即日起到6月30日为止，以优惠价格购买ÄKTAstart主机即可获得相关耗材“红色”大礼包。

800-810-9118

更多关于ÄKTA start 的信息请访问：
www.gelifesciences.com/aktastart

订购信息

| 货号 | 描述 | 数量 | 优惠价 |
|--------------------|--|----|--------|
| 1403Start-tag | 该货号包含以下产品: | | \$8688 |
| | 29-0220-94 ÄKTA start | 1 | |
| | 29-0588-03 Tagged Package, His | 1 | |
| | 29-0588-04 Tagged Package, GST | 1 | |
| | 11-0004-58 HisTrap FF Crude 5 x 1 ml | 1 | |
| | 17-3712-05 HisTrap excel 5 x 1 ml | 1 | |
| | 17-1408-01 HiTrap Desalting 5 x 5 ml | 1 | |
| | 17-6002-33 HiTrap IEX Selection Kit | 1 | |
| 1403Start-antibody | 该货号包含以下产品: | | \$9688 |
| | 29-0220-94 ÄKTA start | 1 | |
| | 29-0588-05 Antibody Package, Protein A | 1 | |
| | 29-0588-06 Antibody Package, Protein G | 1 | |
| | 17-1408-01 HiTrap Desalting 5 x 5 ml | 1 | |
| | 17-6002-33 HiTrap IEX Selection Kit | 1 | |
| | 17-1166-01 HiPrep 16/60 Sephacryl S-200 HR | 1 | |
| | 11-0034-93 HiTrap MabSelect Sure 5*1ml | 1 | |
| 1403Start-Untagged | 该货号包含以下产品: | | \$8688 |
| | 29-0220-94 ÄKTA start | 1 | |
| | 29-0588-07 Untagged Package | 1 | |
| | 17-6002-33 HiTrap IEX Selection Kit | 1 | |
| | 28-9343-88 Copto IEX selection kit | 1 | |

可选模块

| | |
|------------------------------------|------------|
| Frac30 | 29-0230-51 |
| UNICORN start 1.0 control software | 29-0187-51 |



imagination at work